# LABEL:

TURPENTINE LINIMENT 50 ml				
Composition	30 III			
Composition:	DATCH NO	MEC DATE.		
Soft soap :	BATCH NO.	MFG. DATE:		
Camphor :				
Turpentine oil:	MFG. LIC. NO.:	EXPIRY DATE:		
Purified water :	NII OF LICE TO SE	DATINI DITIL		
Category: Counter irritant &				
Rubefacient				
FO	R EXTERNAL USE (	ONLY.		
<b>Use:</b> Liniment of turpentine is	used externally in a pati	ent suffering from arthalgia,		
myalgia, fibrositis and sprain.	_	_		
<b>Storage:</b> Stored in a well closed	d container dark in a co	ol place and temperature not		
exceeding 25°C	a container dans in a co	or place and temperature not		
exceeding 25°C				
MFG BY:	Batch:	Roll No.		

# **Calculation:**

## Experiment - 12

Aim: To	prepare and	submit 50ml	of Tur	pentine	Liniment.
---------	-------------	-------------	--------	---------	-----------

Reference:	 	 

### Requirements:

### **Apparatus Required:**

- Mortar & pestle
- Weigh balance
- Measuring cylinder

Glass rod

Spatula

Beaker

Funnel

### **Chemical Required:**

• Soft soap, camphor, turpentine oil, purified water.

#### Principle:

Liniments are solution or suspensions or emulsion intended for external application. They are generally applied with massage. Liniment of turpentine is used externally in a patient suffering from arthalgia, myalgia, fibrositis and sprain.

### **Turpentine Liniment B.P. Formulation:**

S. No.	Ingredients	Official formula	Required quantities	Uses
1.	Soft soap	9 gm		
2.	Camphor	5 gm		
3.	Turpentine oil	69 gm		
4.	Purified water	1000 ml	50 ml	

#### Procedure:

- ❖ Take the required quantity of soft soap in the mortar and add water in thrice the quantity as soft soap.
- Triturate to make a soapy solution.
- ❖ Take the required quantity of oil of turpentine in a dry measure glass and dissolve camphor in it.
- ❖ Add this solution drop by drop in the mortar and triturate continuously and rapidly till the primary emulsion is formed.
- ❖ Add a small quantity of water and transfer it to the previously calibrated round vertically ribbed and blue and amber colored bottle.
- ❖ Adjust to the required volume by adding water, attach the cork and label it.

## **Marketed Formulation:**









Category: Counter irritant & rubefacient.

#### Use:

- ✓ Liniment of turpentine is used externally in a patient suffering from arthalgia, myalgia, fibrositis, and sprain.
- ✓ It is beneficial in treating pain caused by rheumatoid arthritis and osteoarthritis. It reduces pain and swelling in joints and muscles.
- ✓ It can be used as a muscle relaxant to ease out the pain and stiffness.
- ✓ It provides relaxation in nerve pain and joint pain.

**Storage:** It should be store in a well closed container dark in a cool place.

**Report:** I have prepared the ......ml Turpentine Liniment according to BP and dispense in neat and clean container with proper labeling.

### **Questions Bank**

- **1.** What is Turpentine Liniment?
- **2.** What is the soft soap?
- **3.** What is the use of Turpentine oil in liniment?
- **4.** What is the meaning of Counter irritant?
- **5.** What is the meaning of Rubefacient?
- **6.** What is the Arthalgia disease?
- **7.** Write is the Myalgia disease?
- **8.** What is the Fibrositis?
- 9. What is the Sprain?
- **10.** What is the difference between suspension and emulsion?